Diamond Shield IE Professional

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Mikroströme - EMS

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***Diamond Shield IE
Professional
Wellness Device
Microcurrents - EMS***

Instruction Manual

Firmware Version 2.1

Index

Introduction.....	41
Important Notices.....	43
Intended Usage.....	44
Feature Identification.....	45
Bringing into Service.....	48
Operation.....	49
Menu Navigation.....	50
• Main menu.....	51
• Standard menu.....	52
• Individual menu.....	54
- Define or change a frequency program.....	57
- Working with chipcards.....	59
- Working with a Master chipcard.....	60
• Settings menu.....	64
Technical Data.....	65
For Disposal.....	67
Accessories.....	68
Declaration of Conformity.....	70

Introduction

The Diamond Shield IE Professional Zapper
The Wellness Device with Microcurrents


- wide frequency range of 0.5 Hz to 2.5 MHz.
- intensity range from small voltages (0.1 V) to high voltages (16 V)
- Impulse discharge
- Wobble functionality
- Modulation of the output frequencies
- Constant current control
- Running frequencies from special chip cards (optional accessory)

The Diamond Shield IE Professional allows you to define own frequency programs and store them on Master chipcards (optional accessories). These chipcards can be played in each Diamond Shield device.

The operation is enabled by a 9V battery or a rechargeable battery.

The device can be used to enhance the well-being.

An earthing cable and an enclosed power socket tester are included in the scope of supply. The good effect of the earthing (so called negative ions) occurs only if the power socket is properly earthed. Be sure to check this with the power socket tester.

Warning

The use of the earthing plug (or any other electrical appliances) to an incorrect power socket can be perilous.

Do not use the earthing cable during a thunderstorm.

Important Notices

Safety Instruction

- Use only the supplied grounding cable. Under no circumstances use a different cable for grounding.
- Under any circumstances never insert a plug into an electrical outlet.
- Always insert first the grounding cable to the Diamond Shield and then connect it to the power outlet.
- Never use this device to anyone who has a pacemaker.
- Never use this device in the eye area.
- Keep the device always clean and dry – wiping it with a damp cloth is permitted.

***Using of the device according to the operating instructions is essential.
DermaVit GmbH & Co KG assumes no liability for misuse.***

Note

In case of feeling a tingle or burning on your skin during the usage, please lower the intensity (using the arrow down key). The efficiency of the application is independent from the above. If you are hypersensitive to metal, we will gladly exchange the metal bracelets with cloth bracelets.

Unpacking

Before using this device, please check the completeness of the supplied accessories:

- Connection cable for bracelets or other electrodes (incl. two pairs of connector adapters)
- 2 metal bracelet electrodes
- Grounding cable and power socket tester
- 9V rechargeable battery and battery charger
- Database chipcard

Intended Usage

1. The device can be used to improve the subjective well-being of people who are open to such an approach.
2. It is applicable as EMS-device.
3. With a connected speaker it can serve for acoustic frequency application (music therapy).
4. The preinstalled Diamond Shield program can be used to balance the chinese meridian system.

Storage and transportation

The zapper and its equipment shall be always stored in dry, closed rooms. If the device is transported or stored in extreme temperatures, a period of at least 2 hours should be observed for acclimatization before handling.

The permissible operating temperature ranges from +15°C to +30°C. During storage or transportation, the ambient temperature can be between -20°C and +70°C. If condensation has been formed during transportation or storage, the device must acclimate for about 2 hours and dried by appropriate circulation. Thereafter, the operation of the device is permitted. Nominal data with tolerances are valid after a warm-up time of at least 30 minutes at an ambient temperature of 23°C. A value without tolerances is a guiding value for the device.

The device is intended for use in clean, dry rooms. They cannot be used in extreme dust, damp, dusty air, danger of explosion or chemical vapors.


Protect the Diamond Shield IE Professional from direct sunlight.

Maintenance

The outside of the device should be cleaned regularly with a soft, non-linting cloth. The Diamond Shield IE Professional should never be cleaned with alcohol or other solvents.

Feature Identification

Front Side

- 1. LED
- 2. LCD- Display
- 3. Numeric keypad

- 4. Arrow keys
- 5. START button
- 6. On/Off switch

- 7. Mute switch
- 8. EXIT button
- 9. OK button

Connection Side

1. Grounding connection
2. Chipcard slot
3. Connection for electrodes (ie bracelets)

Back Side

1. Battery compartment

Bringing into service

Insert the battery


1. Open the battery compartment cover on the back of the device.
2. Insert the supplied battery. Please keep in mind that the battery has to be inserted correctly. Meaning the plus has to go to plus and minus to minus. (If however the battery is inserted incorrectly, the device will not be damaged by that.)

We recommend using alkaline batteries or rechargeable batteries.

Connecting the cables

Connect the electrode cable for the bracelet electrodes in the right socket (above LED). Then connect the bracelets.

If you want to work with grounding: Connect the grounding cable to the left socket, then to the outlet.

Operation

Switch On: by pressing briefly the **POWER** button

Switch Off: by pressing the **POWER** button and holding it for about 1 second

Automatic Switch Off: After 12 minutes without user or program activity the device switches off automatically to save battery power.

Charge condition of the battery

The LED indicates the battery charge level. If the device is switched on, the LED flashes. When the battery is full, the flashing frequency is slow. When the battery is almost empty, the flashing frequency is faster.

However the device is fully functional even if the battery is almost empty.

If the battery is completely empty, a message appears, and the device can no longer be used.

Activate / deactivate the loudspeaker

Press the **SOUND** button to switch the audio signal on or off. The current state is then displayed briefly.

Keypad lock


The device has a keypad lock to prevent accidental operation of buttons.

To activate or deactivate the keypad lock briefly press the **POWER**-button and then press **OK**. The device will display the status of the keypad lock.

Menu Navigation

Navigation through the menu structure of the Diamond Shield IE Professional is simple and identical for all functional areas:

The currently selected item is marked with a small black rectangle at the left edge of the display. A scroll bar on the right edge appears, when the current list has more entries than can be displayed. The scroll bar indicates, where you are in the list.

Use the **up-arrow-key** **+** and the **down-arrow-key** **-** to navigate to the desired item.

For switching quickly between the first and the last entry you can jump from the last to the first item with the **down-arrow-key**, and accordingly from the first to the last item with the **up-arrow-key**.

Confirm your selection with the **OK** button. Then a submenu appears, that corresponds to the selected item.

Press the **EXIT** key to return one level up (if you are not already in the main menu).

Important Terms

The structure of the application programs of the standard menu or on chipcards is composed of two levels of elements:

- **Frequency program:** A frequency program is a defined single frequency application with all its parameters such as wave form, wobble, intensity, modulation and duration.
- **Sequence:** A sequence is a succession of several single frequency programs. For example, there is a **sequence** with the name **Diamond Shield**, which includes the **frequency programs DS1, DS2, ... up to DS26**.


Main Menu

Main menu
Chipcard
Standard
Individual
Settings


Shortly after switching on the device, the **Main menu** will appear. Using the **+** and **-** button you can navigate to a menu item and activate the corresponding submenu by pressing the **OK** button. If you want to leave the submenu, you can always return to the **Main menu** by pressing the **EXIT** button.

The Main menu shows the following items:

- **Chipcard:** This item is displayed only if a chip card is inserted. It can be used to display and run or (for Master chipcards) to edit the content of the chipcard (see page 59).
- **Standard:** Under this menu item you can find preinstalled sequences and frequency programs (see page 52).
- **Individual:** This menu allows you to define and edit up to 20 single frequency programs for storing them as a sequence on a Master chipcard. The programs can be executed singularly from this menu (see page 54).
- **Settings:** In this menu various settings (language, display contrast, constant amps, speaker) can be undertaken (see page 64).

Standard Menu


Under the menu item **Standard** you will find the pre-installed frequency applications, with which the Diamond Shield IE Professional Zapper is supplied by default. Most of these applications consist of a sequence of several sub-programs. The first level of the menu shows the available sequences. When a sequence is selected and then the **OK** button is pressed, the frequency programs of the selected sequence are displayed and can be executed singularly.

Executing an entire sequence

Select the desired sequence with the **+** and **-** button. To run the entire sequence press the **START** button.


First the device will test the electrical contact (see chapter **Contact control**, page 53). Then all the sub-programs of the selected sequence are successively executed.

Running a single program of a sequence

Select the desired sequence with the **+** and **-** button. Press the **OK**-button to enter the submenu, where the single frequency programs of this sequence are displayed.

Then select the desired single program and run it by pressing the **START**-button. Previously the device checks the electrical contact (see **Contact control**, page 53).

Contact control

Before a sequence or a single program is executed, the device checks whether the required electrical contact is made. If the contact is not good enough, you will hear a beep every 3 seconds. By pressing the **OK**-button, the frequency program can be executed anyway. However no effect will occur.

If the required contact is made, the display shows for a short time the resistance in kilo-Ohms and then the device automatically begins executing the frequency program.

During the execution of frequency programs

Diam. Shield	DS 1
6:02	5:56
0.6V 0.30mA	0.2V

- In the headline the sequence name and the name of the current single program are shown in turns.
- Below the remaining execution time is counted down (minutes: seconds).
- Below the intensity is displayed.
 - For frequency programs with constant current control, the actual voltage is displayed on the left side. On the right side the chosen intensity is displayed in milliamps (1/1000 amps). (see left figure)
 - For frequency programs without current control, only the voltage in volts is displayed. (see right figure)

Operating options

Pressing the **START**-button will interrupt the program.

Repressing the **START**-button will play the program from the point of interruption.

The intensity can be adjusted using the arrow keys **up** and **down**.

Pressing the **EXIT**-button will quit the current application.

Individual Menu

This menu serves as a working environment. It allows you, to define and edit up to 20 frequency programs, in order to store them as a sequence on a Master chip-card. The stored frequencies are conserved even when the device is switched off or the battery is changed. The frequency programs can be executed singularly by pressing the **START**-button.

You can define the frequency programs with all their parameters using this menu or you can copy programs from a Master- or a database-chipcard and, if necessary, edit them.

In copying on a Master chipcard all contained frequency programs will be included in the sequence. So, for a new composition you have to delete unneeded programs. If you want to retain the defined programs for later compositions, you can use a big Master chipcard quasi as your own database chipcard.

```
Individual
ind. prg. 1
ind. prg. 2
  ---
New freq.
```

The menu **Individual** displays a list of all programs which have been already defined. (Of course for a new device this list will be blank.) The entry **New freq.** allows you to define a new frequency program. see chapter **Define or change a frequency program**, page 55)

```
ind. prg. 1
Edit
Delete
Move up
Move down
```

By selecting a frequency program and pressing the **OK**-button you will enter a submenu, which shows several options for the selected frequency program.

The following actions are available:

- **Edit:** The frequency and all parameters can be edited. See chapter **Define or change a frequency program**, page 55.
- **Delete:** The frequency program is deleted.
- **Move up:** The frequency program is shifted one position up in the list of the menu **Individual**.
- **Move down:** The frequency program is shifted one position down in the list of the menu **Individual**.

Define or change a frequency program

This menu allows you to define a new frequency program or to edit an existing program.

Following entries are available:

<pre>Edit Program name: ind. prg. 1 frequen. 1: 1000.00</pre>	<p>Name: The name of the frequency program, which will be displayed in the menu Individual or on the Master-chip-card. The entry of the name is described in the chapter Name input, page 58.</p> <p>Frequency 1: Fundamental frequency of the program in Hz (the value must range between 0,5Hz and 2.500.000 Hz). A new number is entered using the numeric keys. The previous number is automatically deleted. If only singular digits are to be deleted, this can be done by pressing the arrow key to the left. It is possible but not necessary to specify a decimal place.</p>
<pre>Edit Program name: ind. prg. 1 frequen. 1: 1000.00</pre>	<p>Offset: Offset of the signal</p> <p>The input is done with the arrow keys left and right. You can switch between the options no (no offset), pos (positive offset) and pos + (positive offset + 0.3V).</p>
<pre>Edit Program ind. prg. 1 frequen. 1: 1000.00 offset: post</pre>	<p>Wobble: Enable wobble mode</p> <p>When wobbling, the device varies continuously and cyclically the output frequency between frequency 1 and frequency 2. The input between yes and no is done with the arrow keys left and right.</p> <p>Frequency 2: Second (higher frequency) for wobbling</p> <p>This input is only available if Wobble is enabled. The 2nd frequency must be higher than the 1st frequency. However, it must not be higher than twice the first frequency. The input is done as described at 1st frequency.</p>

```
Edit Program
frequen. 2:
2000.00
wobb. Time:
7.5 sec
```

Wobble Time: Duration of a complete wobble cycle (frequency 1 → frequency 2 → frequency 1). This input is only available if **Wobble** is enabled. The input is done with the arrow keys **left** and **right**. The time must be between 6 and 30 seconds.

```
Edit Program
Const amps:
yes
intensity:
1.00 mA
```

Cons. amps: activation of constant current control
This input is not available when **Const amps** has been disabled in the **Settings**.

In the constant current control the device maintains the current flowing through the body constant. For this purpose, the voltage is increased or decreased according to the instantaneous resistance of the skin and the electrodes. The input is done with the arrow keys **left** and **right**.

```
Edit Program
Const amps:
yes
intensity:
1.00 mA
```

Intensity: Intensity in milliamps (1/1000 amps) or volts, depending on whether the constant current control is enabled or not. The input is done with the arrow keys **left** and **right**.

```
Edit Program
intensity:
1.00 mA
duration:
5:15
```

Duration: Execution time of the program in minutes and seconds. Using the arrow key to the **left** the current time can be deleted. A new time can be entered by using the **numeric keys**. The comma key is used to separate minutes and seconds. If the comma key is not pressed, the number entered will be used in seconds.

For Example: The input of 5,15 results in a time of 5 minutes and 15 seconds.

```
Edit Program
modul.: yes
mod.factor:
4
impul.: no
```

Modul.: enables or disables modulation. The input is done with the arrow keys **left** and **right**.

Mod.-factor: modulation factor between 2 and 254. This command is available only when modulation is enabled. The (multiplication) product of frequency and modulation factor must not exceed 2.500.000 Hz. E.g. at a frequency of 100.000 Hz, a 25-fold modulation is just permissible, a frequency of 500.000 Hz may be modulated only up to 5-fold, a frequency of 9842 Hz can be modulated by 254, thus approaching the modulation threshold (2.499.868 Hz).

A new number is entered using the numeric keys. The previous number is automatically deleted. If only individual digits are to be deleted, this is possible by pressing the arrow key to the **left**.

```
Edit Program
impul.: yes
time act.:
6.00 sec
time pause:
```

Impul.: enables or disables impulse discharge. The input is done with the arrow keys **left** and **right**. If **impulse** is activated, the following commands are available:

Active time: time (in seconds) during which the output frequency is active. The input is done with the arrow keys **left** and **right**.

```
Edit Program
time act.:
6.00 sec
time pause:
1.00 sec
```

Pause time: time (in seconds) when the frequency output is inactive. The input is done with the arrow keys **left** and **right**.

During the defining process the program can be executed by pressing the **START**-button.

Now this program is defined completely.

```
OK: save
EXIT:discard
```

If you press the EXIT button to leave the editing menu, the device asks whether you want **to save the program with OK** or **to discard with EXIT**.

Name input

```
Edit Program
name:
ind. prg.
frequen. 1:
1000.00
```

When saving individual frequency programs or copying a sequence on a Master chipcard, a name can be entered.

The following describes the possibilities:

The name field offers space for 11 letters.

The device already suggests a name. The existing name is automatically cleared when the first letter of a new name is entered.

Should the given name only be completed, then this can be done pushing the arrow key to the **right** before entering the first character. The cursor will jump behind the given name and the automatic overwriting by entering the first letter is suppressed. If single letters should be erased, then this can be done with the arrow key to the **left**.

The input of the text is done with the numeric keys by repeatedly pressing the different keys – similar to the text input with simple mobile phones without T9 predictive text technology, e.g. with the **2**-key: a,b,c,A,B,C,2.


Working with chipcards

As accessory, three different types of chip cards are available:

- Predefined chipcards contain permanently stored sequences and frequency programs. These can be played but neither modified nor copied into the individual memory.
- Master chipcards are delivered empty and can be programmed with user-defined sequences from the individual memory.
- Database-chip cards contain predefined frequency programs, which you can copy to the individual memory for editing and copying from there to a Master Chip Card.

Running a chip card

Regardless of the type of chip card, their sequences and frequency programs can be played as follows:

Insert the chipcard into the chipcard slot of the device.

Select the menu item **Chip card** and press the **OK**-button. A submenu appears, which displays the contents of the chipcard.

The entire sequence(s) and/or single sub-programs can be applied similarly as it is described in the chapter **Standard menu**, page 47.

For running the entire sequence of the selected program, press the **START**-button.

For running single sub-programs of the sequence, enter the submenu of the sequence by pressing the **OK**-button.

The subprograms of this sequence are displayed. Then select the desired subprogram and run it by pressing the **START**-button.


Working with a Master chipcard

Copying all frequency programs from the individual memory to a Master chipcard

The individual frequency programs, which are stored in the individual memory, can be copied to a Master Chip Card. On the Master Chip Card a new sequence will be created and all frequency programs from the individual memory will be stored in that sequence.

Define one or more individual frequency programs (see chapter **Define or change a frequency program**, page 55).


Insert a Master chipcard into the chipcard slot. Select the menu item **Chipcard** from the **Main menu** and press the **OK**-button.

Chipcard
Sequence1

Indiv.->CC
Delete

A submenu will appear displaying the current contents of the chipcard and the additional menu items **Indiv→CC** and **Delete**.

Enter seq.
name:
Sequence2

Select the menu item **Indiv→CC** and press the **OK**-button. A new screen will appear where you can enter the name of the new sequence. Please do this as described in the chapter **Name input**. Then follow the instructions and press the **OK**-button.


All frequency programs from the individual memory are now stored as a new sequence on the Master chipcard.

Deleting the content of a Master chipcard

The entire contents of a Master chipcard can be deleted as follows:

Insert a Master Chip Card into the device.

Select the menu item **Chipcard** from the main menu and press the **OK**-button.

A submenu will appear, displaying the current contents of the Master chipcard and the menu items **Indiv→CC** and **Delete**.

Select the menu item **Delete** and press the **OK**-button.


Confirm the deletion by pressing the **OK**-button.

Depending on the size of the Master Chip Card this process may take up to several minutes.

Rename, move or delete a sequence of a Master chipcard


A sequence stored on the Master chipcard can be renamed, deleted or moved.

Insert a Master chipcard into the chipcard slot. Select the menu item **Chipcard** from the **Main menu** and press the **OK**-button.

A submenu will appear displaying the current content of the chipcard and the additional menu items **Indiv→CC** and **Delete**.

Select the desired **sequence** and press the **OK**-button.

A new menu will appear displaying all frequency programs of the selected sequence and the additional menu items:

- **Seq.→Indiv.**
- **Rename**
- **Move up** (if the selected sequence is not the first one)
- **Move down** (if the selected sequence is not the last one)
- **Delete Seq.**

Select the desired operation and follow the instructions.

Copying a sequence from the Master chipcard to the individual memory

A sequence, which is stored on a Master chipcard, can be copied to the individual memory. In this case, all frequency programs of the individual memory will be overwritten.

The frequency programs of this sequence can then be changed in the individual memory and – if desired – be later copied back to a Master chipcard.

Insert a Master chipcard into the chipcard slot.

Select the menu item **Chipcard** from the **Main menu** and press the **OK**-button.

```
Chipcard
Sequence2
---
Indiv.->CC
Delete
```

A submenu will appear displaying the current contents of the chipcard and the additional menu items **Indiv→CC** and **Delete**.

Select the desired sequence and press the **OK**-button.

```
CC->Seq.2
Seq.->Indiv.
Rename
Move up
Delete Seq.
```

A submenu will appear that displays all the frequency programs of the selected sequence and the additional menu items:

- **Seq.→Indiv.**
- **Rename**
- **Move up** (if the selected sequence is not the first one)
- **Move down** (if the selected sequence is not the last one)
- **Delete Seq.**

```
OK: replace
individual
with
sequence.
EXIT: cancel
```

Select **Seq.→Indiv.** and then press the **OK**-button.

Confirm the prompt by pressing the **OK**-button.

Copying a frequency program from a Master- or a database chipcard to the individual memory

A frequency program can be copied from a Master- or a database chipcard into the individual memory. Here, the individual memory is not overwritten, but the new frequency program is added.

Insert a Master- or a database chipcard into the chipcard slot.

Select the menu item **Chipcard** from the **Main menu** and press the **OK**-button.

```
Chipcard
-----
Sequencel
  ---
Indiv.->CC
Delete
```

A submenu will appear displaying the current contents of the chipcard and the additional menu items **Indiv->CC** and **Delete**.

Select the desired **sequence** and press the **OK**-button.

```
CC->Seq.1
-----
FreqProg1
FreqProg2
FreqProg3
FreqProg4
```


A menu will appear that displays all the frequency programs of the selected sequence and some additional menu items.

Select the desired frequency program and press the **OK**-button.

```
OK: copy
freq. to
individual
EXIT: cancel
```

Confirm the prompt by pressing the **OK**-button.

Settings Menu

This menu is used to configure the device. The settings made here are retained even after switching off the device or changing the battery.

- **Language:** Select from English or German.
The input is done with the arrow keys **left** and **right**.
- **Contrast:** Adjusts the contrast of the LCD display.
The setting between **3** and **20** is done with the arrow keys **left** and **right**.
- **Const amps:** enables or disables the constant current control.
The input is done with the arrow keys **left** and **right**.
If **no** is selected, the intensity of all the frequency programs (including those, that are stored on the chip cards) is usually set on the voltage and does not have a constant current control.
- **Speaker:** the tone pitch of the internal audio signal.
Use the arrow keys **left** and **right** to choose 4 or 5 kHz.
- **HW Rev.:** Displays the hardware version of your device.

Technical Data

General

- Operating temperature: +15° C to +30° C
- Permissible storage temperature: -20° C to + 70° C
- Power supply: 9V alkaline (manganese) battery or a comparable rechargeable battery.

Connections

- 3.5 mm socket (stereo) for grounding plug
- 3.5 mm socket (stereo) for signal outputs
- Chip card slot

Characteristic of the output signals

- Frequency range 0.5 Hz - 2,5 MHz
- Waveform: square wave signal
- Minimum accuracy of the frequencies: 100 ppm
- Output Voltage: 0.1 - 16 Volts (measured from the 0-Volt-line)
- Offset: No offset, positive offset, positive offset + 0.3 Volts
- Modulation up to 254 times the output frequency (but 2,5 MHz must not be exceeded)
- Wobble (maximum up to twice the fundamental frequency)
- Voltage control to keep the output current constant (only at frequencies above 20 Hz)
- Internal resistance: 1 kOhm

For disposal

The Diamond Shield IE Professional device – as well as other electronic devices – must not be disposed with normal household waste. You can bring it to the local separate collection system for electrical appliances or return it to the address below (unfree postage will not be accepted).

Please also note, that the disposal of batteries is prohibited on the normal household waste for environmental reasons. You can bring them to the public collection places or wherever batteries are sold.

You are welcome to return the battery to the following address (unfree postage will not be accepted):

DermaVit GmbH & Co. KG
Unterer Anger 15/Rgb.
80331 Munich
Germany


Accessories

Included in the scope of delivery:

Double cable to connect the electrodes to the Zapper

Two pairs of adapters belong to the electrode cable. The big adapters (10mm) are compatible to all electrodes excepting the adhesive electrodes. These can be connected with the small adapters (4mm). The cable with the red plug is the current carrying line (positive pole).

Metal Bracelet Electrodes

The width of the metal bracelet electrodes can be adjusted. For this purpose the blue element can be opened at the indicated edge.


Database chipcard TOP 121 including the most important frequencies according to Dr. Clark and Dr. Rife

Battery Charger

With the enclosed battery charger you can charge exclusively NiMH (nickel-metal hydride) batteries.

Earthing cable

The insert pins of the earthing cable plug are non-conductive. Only the two lateral earthing contacts are conductive. The cable must first be connected to the zapper, afterwards to a grounded power outlet. The grounding capacity of the socket can be checked by means of the supplied socket tester.

Optional:

Textile Wristband Electrodes

The wristband electrodes' width can be adjusted. The wristbands need to be moistened for use.


Hand Electrodes, stainless steel

The hand electrodes are versatile and easy to clean.

Adhesive electrodes

The adhesive electrodes are available in two sizes, 80 x 45 mm and 35 x 45 mm. They can be placed on the body for specific local application.

Plate, stainless steel

For plate zapping, the plate is interconnected within the current carrying line (the cable with the red plug) by means of an additional cable.

Zappicator

The zappicator is connected to the zapper with the current carrying cable. Then you can place food onto the plate. In the Standard menu you will find the zappicator program. The earthing cable is not needed.

Speaker

The speaker can be used for acoustic frequency application (music therapy).

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© September 2013 DermaVit GmbH & Co. KG
All rights reserved
Any kind of reproduction only with
approval of DermaVit GmbH & Co. KG

DermaVit GmbH & Co. KG
Unterer Anger 15/Rgb.
80331 München
Germany

Tel: +49 (0) 89-18 90 84 6-0
Fax: +49 (0) 89-18 90 84 6-15
E-Mail: info@derma-vit.com
Web: <http://www.derma-vit.com>